

Objective assessment of optic neuritis using chromatic multifocal pupillometry

Zach Weinerman

MS 3, Sackler School of Medicine NY Program

Ygal Rotenstreich Lab

The Maurice and Gabriela Goldschleger Eye Institute,
Sheba Medical Center, Tel-Hashomer

Subjective perimetry



Goldmann




Humphrey

Limitations of subjective perimetry

- Patient cooperation
- More monitoring required
- Cannot distinguish between retina and optic nerve pathologies
- Test-retest variability

Perimetry based on pupillary light reflex to multifocal chromatic stimuli



- ✓ Objective
- ✓ More informative
- ✓ Applicable to various pathologies and patients



Cell Type	Stimulus
Cones	Low-intensity red (624nm)
Rods	Low-intensity blue (485 nm)
ipRGCs	High intensity blue (485 nm)

Objective Multifocal Chromatic Pupillometry (OMCP)

15 parameters




PPC - % pupil contraction

LMCV – Latency of maximal contraction velocity

MCV - Maximal contraction velocity

Exams Measurement Setup

#	X,Y	RGB	Duration	Retest	Data	Change	Baseline
1	-1, 1	0 0 9000	1000ms	No	x	No Data	No Data
2	1, 1	0 0 9000	1000ms	No	x	No Data	No Data
3	1, -1	0 0 9000	1000ms	No	x	No Data	No Data
4	-1, -1	9000 0 0	1000ms	No	x	No Data	No Data
5	-1, -1	9000 0 0	1000ms	No	x	No Data	No Data
6	1, -1	9000 0 0	1000ms	No	x	No Data	No Data



Playback Adjust LED Other

Sequence Control

Enable Sounds

EYE: OD

Name: CHENG, C

MRN: 123

Baseline Pupil Diameter: -



Interval Starting Avg: -

Camera Controls Advanced

Start Camera Stop Camera Save Single Frame

Binary Threshold: 65

--- Measurements Complete
--- Automatic Retests



Warning: Device Hardware Fault Reported

OMCP research

- Previous studies have demonstrated proof of concept for detection of VF defects in patients with:
 - Retinal degeneration
 - Macular dystrophies
 - Glaucoma


- Can OMCP also be effective for detecting optic neuritis?

Optic neuritis

- Unilateral acute inflammation of the optic nerve
- Manifestation of multiple sclerosis
 - Sudden painful vision loss



<http://psych.ucalgary.ca/PACE/VA-Lab/Brian/acquired.htm>



<https://stanfordmedicine25.stanford.edu/the25/fundosopic.html>

Study goal

Utilizing OMCP to assess pupillary responses and identify subtle changes that characterize optic neuritis in the acute phase and during recovery


Study design

- Open, prospective
- 30 acute optic neuritis patients
- 30 MS patients with a history of optic neuritis
- 30 healthy age matched controls

Current enrollment

- 3 acute optic neuritis patients
 - 30 to 38 years old
 - No color vision deficiencies
 - BCVA: 20/30 or better
- 6 healthy age-matched volunteers
- Data collection
 - Within 24 hours of acute ON symptoms
 - After one week of steroid treatment

LMCV - Latency of maximal contraction velocity



Summary

- This preliminary analysis demonstrates that ON patients present with significantly different pupil response kinetics to chromatic light stimuli
- Suggesting that the chromatic multifocal pupillometer may present a novel objective functional test for assessing ON patients

Future directions

- Complete enrollment and analysis of acute ON patients, MS patients and age-matched controls
- Examine the correlation with subjective visual field and optical coherence tomography (OCT) findings
- Determine the sensitivity and specificity of this novel objective visual field testing method for assessing ON disease.



Acknowledgements



Ygal Rotenstreich Lab

- Ifat Sher
- Daniel Ben Ner
- Maya Gurevich
- Helena Solomon
- Irena Kiraj
- Reham Naserat
- Luba Biniaminov
- Lori Gueta

Collaboration:

- Prof. Anat Achiron, Tel Hashomer, Israel

- The Arrow Project
- Sackler School of Medicine NY Program Research Grant